

NEW

Accessory Information

Levigator Stirrup II



7 Reasons to Select Mizuho Stirrups

Interchangeable boots design saves storage space and cost

Quick attachment / release with integrated clamps

Easy and smooth adjustment with gas cylinder and ergonomic handle

Wide range of motion

Lithotomy: +80° to -35°
Abduction: +30° to -10°

**Light-weight design
6.5kg* per unit**

*regular boot, pad and clamp included

Contoured boots with soft and thick pads

Made-in-Japan quality



For more information

For requests and inquiries, please contact our sales company or local sales representative.

Europe | Africa | Middle East

MIZUHO Corporation
EMEA OR Equipment & Device Department
www.mizuho.co.jp/en/mizuho-emea/

Contact us



Asia | Oceania

Mizuho Medical Co.,Ltd.

www.mizuhomedical.co.jp

Contact us



Latin America

Mizuho America, Inc.

www.mizuho.com

Contact us



Regular Stirrup Set

Main Unit + Boot + Pad



	Product Code	Description
R	1 8 - 0 7 0 - 9 1	Levigator Stirrup II Set, Regular, Right Leg
L	1 8 - 0 7 0 - 9 2	Levigator Stirrup II Set, Regular, Left Leg

*Sold separately

Maximize your patient's experience with minimum storage space and investment



- A variety of boot sizes and shapes accommodate different patient types
- Interchangeable boots design saves storage space and cost

Regular

Standard adult size



	Product Code	Description
R	1 8 - 0 7 0 - 2 1	Boot for Levigator Stirrup II, Regular, Right Leg
L	1 8 - 0 7 0 - 2 2	Boot for Levigator Stirrup II, Regular, Left Leg
R	1 8 - 0 7 0 - 4 1	Boot Pad for Levigator Stirrup II, Regular, Right Leg
L	1 8 - 0 7 0 - 4 2	Boot Pad for Levigator Stirrup II, Regular, Left Leg

*Sold separately

Fin Type

Lateral knee support further reduces risk of nerve damage



	Product Code	Description
R	1 8 - 0 7 0 - 2 3	Boot for Levigator Stirrup II, Fin Type, Right Leg
L	1 8 - 0 7 0 - 2 4	Boot for Levigator Stirrup II, Fin Type, Left Leg
R	1 8 - 0 7 0 - 4 3	Boot Pad for Levigator Stirrup II, Fin Type, Right Leg
L	1 8 - 0 7 0 - 4 4	Boot Pad for Levigator Stirrup II, Fin Type, Left Leg

*Sold separately

Wide Type

For large, overweight and obese patients



	Product Code	Description
R	1 8 - 0 7 0 - 2 7	Boot for Levigator Stirrup II, Wide, Right Leg
L	1 8 - 0 7 0 - 2 8	Boot for Levigator Stirrup II, Wide, Left Leg
R	1 8 - 0 7 0 - 4 7	Boot Pad for Levigator Stirrup II, Wide, Right Leg
L	1 8 - 0 7 0 - 4 8	Boot Pad for Levigator Stirrup II, Wide, Left Leg

*Sold separately

Main Unit



	Product Code	Description
R	1 8 - 0 7 0 - 9 3	Levigator Stirrup II Main Unit, Right Leg
L	1 8 - 0 7 0 - 9 4	Levigator Stirrup II Main Unit, Left Leg

*Sold separately

Levigator Cart



Product Code	Description
1 8 - 0 7 0 - 6 1	Levigator Stirrups Cart